

ICE FLAKER



Ice Flaker

FEATURES

- Fully Automatic microprocessor control
- Compact Design and low maintenance
- Stainless Steel construction
- Continuous Ice Flakes Output
- CFC free compressor
- Noiseless operation
- Overload protection
- Low water level Indication

APPLICATION :

Research Institute, Hospitals, Laboratories
Biotechnology & Universities,

SPECIFICATIONS

MODEL	LC-57-A	LC-57-B	LC-57-C	
Capacity of Ice flake (kg per 24 hrs.)	30 kg	50 kg	100 kg	
Storage Capacity (kg)	10 kg	15 kg	25 kg	
Storage bin	ABS Plastic			
Cooling Method	Air			
Ambient temperature	Below 30 °C for optimum working			
Construction	Stainless Steel			
Noise level	≤52 dB			
Water Supply	3 / 4 inch BSP			
Refrigerant	R134A			
Input Power	180W	280W	420W	
Dimension (mm)	300X493X547	380X543X722	548X611X883	
Power supply	220V / 50Hz			
Gross Weight	30 kg	36kg	45kg	64kg

Note : High Capacity models are also available - 150kg, 200kg, & 300kg.

* Specifications subject to change without notice.