



MANUFACTURERS OF:
Industrial Oven, Furnaces, Cooling
Equipments, Scientific, Surgical, &
Laboratory Equipments

LAB-CARE
INSTRUMENTS & SERVICES

Blood Bank Tube Sealer Unit

LAB CARE Blood Bag Tube Sealer is a compact equipment to seal the Blood Bag pilot tube by radio frequency sealing system. Sealing can be achieved in less than 1.5 seconds per tube by simple operation. Sealed tube can be easily separated by pulling the sealing section by both sides

Technical Specification	
Model	LC-62
Oscillation Frequency	40.68 MHz.
Radio Frequency Output	15 W
Sealing time (Approx.)	1.5 Sec
Input Voltage	230V AC
Frequency	50 Hz.
Power Consumption	100V A
Weight Approx.	5.5 Kgs.



*N.B. : Due to continuous R&D model / specifications may change.
Picture shown with optional accessories.*